



Website: <http://performhealth.com.au/Product/info/OPP-2264-12S>

* Designed with straps crossing at the back and fastened in front. * Comfortable to wear and provides excellent support. * Four pliable support splints at the back for extra stability. * Breathable elastic and soft pile bonded to foam. * Available in two back heights: 9" (22.8cm) and 12" (30.4cm) How to wear: * Wrap support around the waist and be sure the four stays are on each side of spine. * Adjust to desired compression and secure the hook and loop strap to the front of support. * Pull the elastic reinforcement straps to tighten, then secure to the front of support. Application: * Low back pain * Muscle strain and ligament sprain * Osteoporosis * Spondylosis

Specifications

Google Category	Health & Beauty > Health Care > Supports & Braces
Product Detail Group	Oppo 2264 Sacro Lumbar Support
Product List Group	Oppo 2264 Sacro Lumbar Support
Height	12"
GOOGLECAT	Health & Beauty>Health Care>Supports & Braces
EBAYCAT	Health & Beauty>Medical, Mobility & Disability>Orthopaedics & Supports>Orthotics, Braces & Sleeves
eBay_SKU	OPP-2264
eBay_Title	Oppo 2264 Sacro Lumbar Support
Size	S